





INSULCAP THERMAL PALLET COVERS

Thermal barrier systems enclose the palletised product to provide an protected thermal environment, that passively manages product temperature during storage and transport. An InsulCap® cover is constructed from a range of thermal barrier materials; innovations designed by Wilpak Group International.

WHY CHOOSE INSULCAP THERMAL COVERS?


- A range of custom InsulCap® pallet covers; InsulPrime, InsulGold, InsulPlatinum Plus
- Our InsulPlatinum® Plus barrier is a heavy duty large air cylinder material which offers maximum protection for temperature
- sensitive goods, while our InsulPrime® barrier is a generalpurpose material at a affordable price.
- Exclusively Distributed Globally by Peli BioThermal

EXCELLENT THERMAL PROTECTION

TOUGH THERMAL BARRIER MATERIAL

LIGHTWEIGHT PALLET COVER SYSTEM

CUSTOM CUT PROFILES AVAILABLE IN 150MM INCREMENTS





INSULCAP THERMAL BARRIER MATERIAL PROPERTIES

	Reflective Barrier	Air Cy	linder	Barrier T	Thermal Rating	
		mm		mm	inches	
InsulPlatinum Plus	2 (foil external/internal)	25+	0.98	12	0.47	5 Maximum
PloBlusul	2 (foil external)	25+	0.98	10	0.39	4 High
InsulPrime	1 (foil external)	10	0.39	5	0.19	1 Affordable

THERMAL BARRIER MATERIAL

REQUIRED TEMPERATURE	EXTERNAL EXPOSURE TEMPERATURE	EXPOSURE DURATION							
RANGE OF GOODS		Hrs	1	2	4	6	8	8+	
Lower Upper Limit Limit		Mins	60	120	240	360	480	480+	
2°C 8°C	less than -10°C							Performance Varies Based on Application	
20 80	-10°C to 0°C								
	0°C to 10°C		INSL	ادیال					
1966 & TEMPERATURE SOMETIVE +2°C to +8°C	10°C to 20°C								
+2*C to +8*C	20°C to 30°C		175		_ATIN	JUM.			
Chilled	30°C to 40°C					PLUS		P B;	
	200 0 000								
15°C 25°C	less than -10°C							se	
700000	-10°C to 0°C			ULP	LATIN	IUM.		/aric	
	0°C to 10°C					PLUS		nce \	
TIME & TEMPERATURE SENSITIVE +15°C to +25°C	10°C to 20°C			-12/5				rmai on /	
710 0 10 725 0	20°C to 30°C							Performance Varies Based on Application	
Ambient	30°C to 40°C							P B:	



FOR MORE INFORMATION ABOUT PELI BIOTHERMAL PRODUCTS **CLICK HERE**

FOR MORE INFORMATION ABOUT PELI BIOTHERMAL PRODUCTS, PLEASE CONTACT:

EMEA OFFICE - INFOEMEA@PELIBIOTHERMAL.COM T: +44 (0)1525 243 770

AMERICAS OFFICE - INFOAMERICAS@PELIBIOTHERMAL.COM T: +1(763) 412-4800 TOLL FREE US: T: +(877) 537-9800

ASIA OFFICE - INFOASIA@PELIBIOTHERMAL.COM T: +65 6681 0095

PELIBIOTHERMAL.COM

